

# Ensuring Quality and Brand Consistency Through Color Measurement of Potato Chips Using a HunterLab Aeros Spectrophotometer

## I. The Critical Role of Color in Potato Chip Quality and Consumer Acceptance

The visual appeal of food products significantly influences consumer perception, and among the various attributes, color stands out as the primary indicator of quality and freshness for potato chips.<sup>1</sup> Research consistently demonstrates that consumers make their initial assessment of potato chips based on their color, directly linking this visual characteristic to their expectations of flavor and overall quality.<sup>3</sup> A desirable light to light golden-yellow hue, often described as golden-brown, is typically associated with properly cooked and fresh potato chips, signaling an appealing taste and texture.<sup>1</sup> Conversely, if potato chips appear too pale, it might suggest undercooking, while an excessively dark color, such as brown or burnt, often leads to negative perceptions and potential rejection by consumers who may associate it with undesirable flavors or lower quality.<sup>1</sup> This immediate impact of visual appeal underscores the critical necessity for potato chip manufacturers to implement precise color control measures throughout their production processes.

Beyond the initial impression, consistent color plays a vital role in shaping consumer preference and fostering brand loyalty.<sup>3</sup> Consumers tend to favor brands that exhibit uniformity in color across different batches of potato chips, suggesting that consistency is a key factor in building and maintaining brand preference.<sup>3</sup> Inconsistent color can lead consumers to perceive quality issues, potentially prompting them to explore alternative brands that offer a more reliable visual experience.<sup>3</sup> Therefore, maintaining a stable and recognizable color profile is essential for reinforcing brand recognition and ensuring that customers feel confident in the quality of the product with every purchase.<sup>8</sup> Color acts as a visual signature of the brand, and variations can erode consumer trust and ultimately impact their loyalty.

Accurate color measurement is also indispensable during the stages of product development and quality assurance.<sup>3</sup> When manufacturers introduce new flavors or modify production techniques, understanding how these changes affect the final color of the potato chips is crucial. Objective color data obtained through precise measurement allows for informed decisions regarding ingredients and processes.<sup>3</sup> Furthermore, monitoring color at each stage of production serves as a proactive quality assurance measure. By tracking color, manufacturers can identify when adjustments are necessary to prevent waste and ensure that only high-quality

products reach the market.<sup>3</sup> This objective data also helps manufacturers adhere to established color standards, which may be internal specifications or driven by regulatory requirements or customer expectations.<sup>3</sup> Color measurement, therefore, is not merely a final inspection step but an integral component of the entire product lifecycle, from initial design to continuous quality control.

While a single set of universal regulations specifically governing potato chip color might not be prevalent, various industry standards and guidelines emphasize the importance of this attribute.<sup>4</sup> For instance, the United States Department of Agriculture (USDA) provides standards for grades of potatoes intended for chipping, which include optional tests for fry color utilizing the Agtron scale.<sup>10</sup> These standards outline acceptable color ranges and recommend specific methods for measurement, highlighting the industry's broad recognition of color as a significant quality parameter. Adherence to these recognized industry standards, even when optional, demonstrates a manufacturer's commitment to producing high-quality potato chips and can contribute to greater consumer confidence and facilitate trade within the industry.

## **II. Objective Color Measurement: The Science of Spectrophotometry and the HunterLab Aeros**

Spectrophotometry provides a scientific foundation for the precise and objective measurement of color in various materials, including food products like potato chips.<sup>9</sup> Spectrophotometers function by quantifying the interaction of light with a sample, specifically by measuring the amount of light that is either absorbed or reflected across a defined spectrum of wavelengths.<sup>9</sup> This analysis yields objective, quantifiable data that characterizes the sample's color properties.<sup>9</sup> Unlike subjective assessments based on human vision, spectrophotometry eliminates the inherent variability associated with individual perception and provides measurements that are both consistent and repeatable.<sup>3</sup> This technology offers a reliable method for color analysis, moving beyond the inherent limitations of human sensory evaluation.

The advantages of employing spectrophotometry over subjective visual assessment in the context of potato chip quality control are numerous and significant.<sup>8</sup> Visual color analysis is susceptible to inaccuracies arising from individual differences in how people perceive color, variations in lighting conditions, and other subjective factors such as observer fatigue.<sup>8</sup> In contrast, spectrophotometers furnish precise numerical data that can be used to establish and monitor color standards with a much higher degree of accuracy.<sup>3</sup> Furthermore, instrumental analysis enables the detection of subtle color variations that might not be perceptible to the human eye, which is

particularly important for maintaining stringent quality control in potato chip production.<sup>3</sup> The adoption of spectrophotometry represents a transition towards a data-driven approach to quality control, resulting in more dependable and consistent product outcomes.

The HunterLab Aeros spectrophotometer is specifically engineered to address the unique challenges associated with color measurement of irregular, non-flat samples, making it particularly well-suited for potato chip analysis.<sup>3</sup> A key feature of the Aeros is its non-contact measurement methodology, which simplifies the process of sample presentation and eliminates the potential for contamination of either the instrument's sensor or the sample itself.<sup>3</sup> Additionally, the Aeros provides standardized data that can be directly compared to widely recognized color indices and scales commonly used within the potato chip industry.<sup>3</sup> This tailored design and functionality make the HunterLab Aeros a valuable tool for ensuring color quality and consistency in potato chip manufacturing.

Several key features and benefits of the HunterLab Aeros make it an exceptional choice for potato chip analysis.<sup>14</sup> Its non-contact measurement capability simplifies sample handling, reduces the risk of contamination, and minimizes the time required for cleanup.<sup>14</sup> The automatic height positioning feature ensures that the sensor maintains an optimal distance from the sample, regardless of its thickness or shape, thereby enhancing the accuracy of measurements.<sup>14</sup> The rotating sample platform provides what is described as the world's largest measurement area for a benchtop spectrophotometer, allowing for the capture of color data from a representative portion of the typically non-uniform potato chip sample.<sup>14</sup> This system can acquire multiple color and spectral measurements across a substantial area in a matter of seconds, significantly improving efficiency.<sup>14</sup> The embedded EasyMatch® Essentials software offers a user-friendly interface and includes essential color indices and metrics such as CIE L\*a\*b\* and Baking Contrast Units (BCU).<sup>14</sup> Furthermore, the Aeros is designed for seamless integration with various quality control methodologies, including Statistical Process Control (SPC) systems, providing versatility in its application within a manufacturing environment.<sup>3</sup> These combined features make the Aeros a comprehensive solution for objective and efficient color measurement of potato chips, directly addressing the practical needs of the industry.

### **III. Standard Operating Procedures for Accurate Color Measurement Using the HunterLab Aeros**

To ensure the reliability and consistency of color measurements obtained using the HunterLab Aeros spectrophotometer for potato chips, a well-defined standard

operating procedure (SOP) is essential. This SOP should encompass instrument setup, calibration, standardization, sample preparation, and data acquisition.

The initial step involves the proper setup of the HunterLab Aeros.<sup>17</sup> This includes unpacking the instrument and placing it on a stable, level surface, followed by connecting the power cable and turning on the power switch located on the back of the instrument base. Once powered on, it is crucial to allow the instrument to warm up for the recommended duration, typically around 30 minutes, to ensure that it functions correctly and provides accurate measurements.<sup>18</sup> After the warm-up period, the EasyMatch® Essentials software should be initialized, and the connection between the software and the instrument should be verified.<sup>17</sup> The status bar within the software interface usually indicates whether the instrument is successfully connected.

The next critical step is standardization, which calibrates the instrument using certified reference standards to ensure the accuracy and traceability of all subsequent measurements.<sup>17</sup> Standardization can be initiated either through the Workspace menu within the software or by clicking on the standardization status area typically found on the status bar. The software will prompt the user to remove any samples from under the sensor. Following the prompts, the sensor will move to its highest position, and the user will be instructed to attach the provided standardization box to the sensor. The black glass reference standard is then attached to the standardization box, and the "READ" button is pressed. After reading the black glass, it is removed, and the calibrated white tile is attached in its place, followed by another press of the "READ" button. Once this process is complete, the calibration box is removed. If the option to read the green tile during standardization is enabled in the software preferences, the user will also be prompted to measure the green diagnostic check tile to further verify the instrument's calibration.<sup>17</sup> Upon successful completion, the status bar will typically update to indicate that the instrument is "STANDARDIZED".

Proper sample preparation and presentation are equally important for obtaining representative and accurate color measurements, especially with non-uniform products like potato chips.<sup>17</sup> A representative sample of potato chips should be collected from the production line or the batch being tested, ensuring that the sample size is sufficient to adequately cover the measurement area of the Aeros.<sup>17</sup> The collected potato chips should then be placed in a clean sample dish, which may be a standard petri dish or a specialized dish provided with the Aeros.<sup>14</sup> It is important to distribute the chips evenly within the dish to create a consistent surface for measurement, avoiding significant overlapping or large gaps as much as possible.<sup>11</sup> A minimum depth of 1 inch is generally recommended to ensure the bottom of the sample tray is fully covered.<sup>11</sup> The prepared sample dish is then placed in the center of

the Aeros' rotating sample platform.<sup>11</sup>

The operation of the Aeros software for data acquisition involves a few straightforward steps.<sup>11</sup> It is advisable to first create a new workspace within the EasyMatch® Essentials software specifically designated for potato chip measurements. This allows for the customization and saving of instrument settings and tolerance limits relevant to this particular product.<sup>17</sup> Within this workspace, the read options should be configured. This includes setting the measurement time, with 5 seconds being the default and often suitable for potato chips, and ensuring that the turntable is enabled ("ON") to facilitate automatic sample rotation during measurement.<sup>17</sup> The Auto Height positioning feature should also be enabled, as this allows the sensor head to automatically adjust its height to the optimal distance for each sample.<sup>17</sup> Once the sample is in place and the software is configured, the measurement process is initiated by pressing the "READ" button on the software interface.<sup>11</sup> The Aeros sensor will then automatically scan across the rotating sample, capturing multiple color measurements.<sup>11</sup> The resulting color measurement data will be displayed on the screen, typically showing the values for the pre-configured color scales, such as CIE L\*a\*b\*.<sup>14</sup> This data should be recorded for subsequent analysis and quality control documentation. The software also provides the option to name each sample, which aids in organization and identification.<sup>17</sup> For enhanced accuracy, it is often recommended to take multiple readings of the same sample, slightly rotating the sample dish between each measurement to account for potential variations within the sample itself.<sup>19</sup> Finally, the measurement data is typically saved within a job file associated with the specific workspace, and the software offers options for recalling, importing, exporting, emailing, or deleting data as needed for quality management purposes.<sup>17</sup>

### **Must-have Table: Key Steps in SOP for HunterLab Aeros Measurement of Potato Chips**

<b>Step</b>	<b>Action</b>	<b>Details</b>
1	Instrument Setup	Place on stable surface, connect power, turn on.

2	Warm-up	Allow instrument to warm up for ~30 minutes.
3	Software Initialization	Launch EasyMatch® Essentials, ensure instrument connection.
4	Standardization	Use black glass and white tile reference standards as prompted.
5	Workspace Creation	Create a dedicated workspace for potato chips.
6	Configure Read Options	Set measurement time (e.g., 5 seconds), enable turntable and Auto Height.
7	Sample Preparation	Collect representative sample, distribute evenly in sample dish.
8	Sample Placement	Place dish on the center of the rotating platform.
9	Initiate Measurement	Press "READ" in the software.
10	Data Recording	Note the displayed color values (e.g., L*a*b*).
11	Multiple Readings (Optional)	Rotate sample and repeat measurement for increased accuracy.
12	Data Management	Save and manage data within the software.

#### **IV. Decoding Color: Relevant Scales and Indices for Potato Chip Assessment**

The accurate assessment of potato chip color relies on the use of appropriate color scales and indices, which provide a standardized way to quantify and communicate color attributes. Several color spaces and indices are particularly relevant in the

context of potato chip quality control, including CIE L\*a\*b\*, Hunter L, a, b, the Agron scale, and Baking Contrast Units (BCU).

The CIE L\*a\*b\* color space, defined by the International Commission on Illumination (CIE), is a widely adopted system for describing colors based on three perceptual attributes.<sup>20</sup> These attributes are lightness (L\*), redness-greenness (a\*), and yellowness-blueness (b\*).<sup>21</sup> The L\* value ranges from 0 (representing perfect black) to 100 (representing perfect white), with higher L\* values in potato chips generally indicating a lighter, more desirable color.<sup>21</sup> The a\* value indicates the degree of redness or greenness; positive values signify redness, while negative values signify greenness.<sup>21</sup> For potato chips, a slightly positive a\* value contributes to the characteristic golden hue, whereas negative values might suggest under-frying, and excessively high positive values could indicate burning or over-processing.<sup>23</sup> The b\* value represents yellowness or blueness, with positive values indicating yellowness and negative values indicating blueness.<sup>21</sup> A positive b\* value is crucial for achieving the desired golden-yellow color that consumers typically associate with high-quality potato chips.<sup>5</sup> One of the key advantages of the CIE L\*a\*b\* color space is its perceptual uniformity, meaning that a numerical change in the color space corresponds to a relatively consistent perceived change in color by the human eye.<sup>20</sup> The HunterLab Aeros software suite typically includes CIE L\*a\*b\* as a standard option for color measurement and data reporting.<sup>14</sup> Understanding the L\*a\*b\* values obtained from the Aeros allows for a precise and objective characterization of potato chip color, facilitating the establishment of specific color targets and tolerance limits for quality control.

The Hunter L, a, b scale is another color measurement system, developed by HunterLab, which has become a foundational element in the field of spectrophotometry.<sup>24</sup> Similar in concept to CIE L\*a\*b\*, it also employs three values to define color: L for lightness (ranging from 0 for black to 100 for white), a for redness-greenness (with positive values indicating red and negative values indicating green), and b for yellowness-blueness (where positive values represent yellow and negative values represent blue).<sup>25</sup> While the underlying principles are comparable to CIE L\*a\*b\*, subtle differences may exist in the mathematical transformations used to derive these values. Given that the HunterLab Aeros is the instrument in question, familiarity with the Hunter L, a, b scale is directly pertinent to interpreting the color data it produces.<sup>14</sup> Although both CIE L\*a\*b\* and Hunter L, a, b aim to represent color in a manner that aligns with human visual perception, it is important to recognize that the numerical values for the same sample might not be identical across the two systems. Therefore, for internal consistency in quality control, it is advisable to

primarily utilize one scale or the other.

The Agtron scale is a colorimetric system specifically designed for use in the food industry, particularly for evaluating the color of fried products such as potato chips.<sup>26</sup> This scale provides a single numerical index that correlates with the perceived degree of browning, with higher numbers on the scale indicating lighter colors.<sup>26</sup> The USDA standards for potatoes intended for chipping reference the Agtron scale and provide specific target index ranges for different color designations, ranging from very light to dark.<sup>10</sup> Agtron analyzers are widely employed in the potato chip manufacturing industry for routine quality control purposes.<sup>27</sup> While the HunterLab Aeros might not directly output Agtron values, it is often possible to establish a correlation or conversion between the L\*a\*b\* or Hunter L, a, b values it measures and the corresponding Agtron scale values.<sup>10</sup> This is particularly useful as the Agtron scale is an industry-specific standard that offers a practical way to communicate and control potato chip color based on established and widely understood categories.

Baking Contrast Units (BCU) represent another color index utilized in the baking industry to assess the extent of browning or contrast in baked goods, which includes potato chips.<sup>14</sup> The HunterLab Aeros software typically includes BCU as one of the color indices available for measurement and analysis.<sup>14</sup> This index is particularly valuable for monitoring the progression of the Maillard reaction and the overall browning process that occurs during the frying of potato chips. By providing a specific metric for quantifying the browning, BCU allows for enhanced control over the frying process to achieve the desired product characteristics in terms of color, flavor, and texture.

## **V. Establishing Color Consistency: Recommended Ranges and Target Values for Potato Chips**

Achieving color consistency in potato chip production necessitates the establishment of clear and measurable color ranges and target values. These standards should be informed by consumer preferences, brand specifications, and relevant industry guidelines.

Consumer studies consistently indicate a preference for potato chips that exhibit a light to light golden-yellow color.<sup>7</sup> A golden-brown appearance is frequently cited as the most desirable color attribute for high-quality potato chips.<sup>1</sup> However, each individual brand may have its own specific target color range that serves to differentiate its product in the marketplace and align with the expectations of its particular consumer base. These color targets can be defined using a variety of color

scales, including the CIE L\*a\*b\* system, the Hunter L, a, b scale, and the industry-specific Agtron scale. It is therefore crucial that manufacturers base their color targets on a comprehensive understanding of both general consumer preferences for potato chip color and the unique visual identity and product positioning of their own brand.

Existing industry standards and reference charts can provide a valuable starting point for defining these color targets. For example, the USDA standards for potatoes for chipping provide Agtron index ranges that correspond to different fry color designations.<sup>10</sup> These designations range from "1," which indicates an Agtron index of 65 or higher and represents a very light color, to "5," which corresponds to an index of 25 to 34, signifying a darker color. While instrumental methods like spectrophotometry offer greater objectivity, visual reference charts, such as the 10-point color chart developed for use in India, can also be utilized as guides for assessing chip color.<sup>7</sup> Additionally, the European Potato Processors Association has identified a light golden color as generally acceptable for fried potato products.<sup>4</sup> These existing standards and reference materials offer benchmarks based on industry best practices and insights into consumer expectations, which can be adapted and refined to meet specific brand requirements.

It is also important to consider that acceptable color ranges may vary depending on the specific style of potato chips being produced.<sup>30</sup> For instance, kettle-cooked potato chips, which are often fried at lower temperatures for a longer duration, may naturally exhibit a slightly darker color compared to traditional thin-cut varieties. Similarly, the addition of seasonings to create flavored potato chips can significantly influence the overall color appearance of the final product. Therefore, when establishing color targets, manufacturers must take into account the inherent characteristics and consumer expectations associated with each particular style of potato chip in their product line. A single, uniform color standard may not be appropriate across all styles.

**Must-have Table: Recommended Color Ranges for Various Potato Chip Styles (Example)**

Potato Chip	Desired Color Descripti	L* Range (Example	a* Range (Example	b* Range (Example	Agtron Index Range	BCU Range (Example
-------------	-------------------------	-------------------	-------------------	-------------------	--------------------	--------------------

Style	Color	L*	a*	b*	(USDA)	Range
Traditional Thin-Cut	Light Golden Yellow	65-75	-2 to 2	30-40	1-2 (65+)	60-70
Kettle-Cooked	Golden Brown	55-65	2-5	25-35	2-3 (55-64, 45-54)	50-60
Lightly Salted	Pale Golden	70-80	-3 to 1	35-45	1 (65+)	65-75
Barbecue Flavored	Medium Brown	50-60	5-10	20-30	3-4 (45-54, 35-44)	45-55

*Note: These are example ranges and actual values should be determined based on specific brand requirements and consumer testing.*

**VI. The Influence of Processing Parameters on Potato Chip Color**

The final color of potato chips is a direct result of the complex interplay of several processing parameters, with frying time and temperature being among the most critical. The characteristic color of potato chips develops primarily through the Maillard reaction, a non-enzymatic browning process that occurs when reducing sugars and amino acids in the potato slices react under the high heat of frying.<sup>1</sup> Generally, as frying time is extended and the frying temperature is increased, the potato chips will become darker in color due to the continued progression of the Maillard reaction.<sup>32</sup> Interestingly, some studies suggest that frying at lower temperatures for a longer duration can sometimes lead to higher L\* values, indicating a lighter color, compared to frying at higher temperatures for a shorter time.<sup>33</sup> However, it is important to note that excessive browning can result in the development of undesirable burnt flavors and an unacceptable dark color in the finished product.<sup>1</sup> Therefore, precise control over both frying time and temperature is paramount for achieving the desired color consistency in potato chip production.

Beyond the immediate frying process, the characteristics of the raw potato itself, including the variety, its maturity at harvest, and the conditions under which it was stored, exert a significant influence on the final chip color.<sup>1</sup> Different potato varieties inherently contain varying levels of reducing sugars, which are the primary reactants

in the Maillard browning process. Potatoes with higher levels of reducing sugars will tend to produce darker chips when fried.<sup>1</sup> The maturity of the potato crop at the time of harvest can also affect its sugar content and, consequently, the resulting chip color.<sup>1</sup> Furthermore, storage conditions, particularly temperature, can have a substantial impact on the conversion of starch to sugars within the potato, a phenomenon known as cold-induced sweetening. Potatoes stored at low temperatures tend to accumulate higher levels of reducing sugars, which can lead to the production of darker chips during frying.<sup>4</sup> To mitigate this effect in cold-stored potatoes, a process called reconditioning, which involves storing the potatoes at higher temperatures for a period before processing, can be employed to help reduce the sugar content.<sup>4</sup> Thus, the selection of appropriate potato varieties and careful management of storage conditions are crucial for minimizing color inconsistencies in the final potato chip product.

Pre-treatment processes applied to the potato slices before frying can also play a role in influencing their final color. Washing the potato slices is an essential step to remove surface dirt and debris, which can indirectly affect the consistency of the frying process and the resulting color. Blanching, which involves briefly immersing the potato slices in hot water, can serve several purposes that impact color. It can inactivate enzymes present in the potato that might otherwise cause discoloration, and it can also leach out some of the reducing sugars from the slices. The removal of these sugars can potentially lead to the production of lighter-colored chips after frying.<sup>34</sup> The effectiveness of blanching in controlling color is dependent on the specific temperature and duration of the treatment. Therefore, careful optimization of these pre-treatment steps can contribute to greater color consistency in the final potato chip product by influencing the composition of the potato slices before they enter the fryer.

## **VII. Implementing Strategies for Maintaining Color Uniformity in Potato Chip Production**

To consistently produce potato chips that meet the desired color standards and brand specifications, manufacturers can implement several key strategies, with the data obtained from the HunterLab Aeros spectrophotometer playing a central role in this effort.

One effective strategy involves integrating the HunterLab Aeros directly into the production line to enable real-time or frequent color monitoring of the potato chips as they are being produced.<sup>3</sup> By providing continuous data on the color of the chips, this approach allows operators to closely track color trends and make immediate

adjustments to critical processing parameters, such as frying time and temperature, should any deviations from the established target range be detected.<sup>40</sup> The software associated with the Aeros can be configured with pre-defined color standards and acceptable tolerance limits. If the measured color values fall outside these limits, the system can generate alerts, prompting operator intervention to correct the process before an entire batch of off-specification product is produced.<sup>40</sup> This real-time color monitoring capability facilitates a proactive approach to quality control, significantly minimizing the production of out-of-spec batches and reducing potential waste.

The valuable color measurement data generated by the HunterLab Aeros can also be integrated into Statistical Process Control (SPC) systems.<sup>3</sup> By incorporating this data into SPC charts and analyses, manufacturers can effectively track the stability of their production process over time and identify potential sources of variation that may be contributing to color inconsistencies.<sup>3</sup> Analyzing trends in color measurements can help predict when adjustments to the process might be necessary to proactively maintain color uniformity. Furthermore, SPC can be utilized to establish statistically derived control limits for key color parameters. These limits provide an objective and data-driven basis for making quality control decisions and ensuring that the process remains within acceptable boundaries. The integration of color data into SPC offers a powerful tool for understanding and managing process variability, ultimately leading to more consistent and high-quality potato chip production.

Establishing clear and well-defined tolerance limits for color measurements is another crucial step in maintaining color uniformity. Based on the established color targets and the specific brand specifications for the potato chips, manufacturers should define acceptable upper and lower limits for the key color parameters, such as  $L^*$ ,  $a^*$ ,  $b^*$ , or Agtron values. Alongside these tolerance limits, it is essential to develop clear and actionable protocols outlining the corrective measures to be taken whenever color measurements fall outside of these defined boundaries. These protocols might include adjusting frying time or temperature, conducting checks on the quality of the incoming potatoes, or investigating other potential factors in the production process that could be contributing to the color variation. Having these pre-defined tolerance limits and corrective action plans in place ensures a systematic and consistent approach to managing color deviations, thereby minimizing subjective judgments and promoting a proactive strategy for upholding product quality.

Finally, achieving consistent color in potato chip production requires a comprehensive approach that addresses all stages of the process, starting with the sourcing of raw materials and extending through to the final frying stage. It is vital to establish relationships with reliable potato suppliers who can consistently provide potatoes of

the desired variety and with stable, predictable sugar content.<sup>1</sup> Optimizing the pre-treatment processes, including washing, slicing, and blanching, is also critical to minimize variability in the potato slices before they are fried. This ensures a more uniform starting material for the frying process. Furthermore, regular calibration and diligent maintenance of the frying equipment are essential to guarantee consistent oil temperature and frying time across all production batches. By focusing on the consistency of the raw materials and the tight control of processing parameters, manufacturers can significantly reduce the likelihood of unwanted color variations in their potato chips and ensure a final product that consistently meets consumer expectations and brand standards.

### **VIII. Real-World Applications: Case Studies of HunterLab Spectrophotometers in Snack Food Quality Control**

The effectiveness of HunterLab spectrophotometers in ensuring color quality and consistency is well-recognized within the snack food industry. HunterLab themselves highlight that many of the world's leading potato chip manufacturers rely on their advanced color measurement solutions for maintaining stringent quality control.<sup>8</sup> The HunterLab Aeros, with its features specifically designed for the challenges of measuring irregular samples like potato chips, is often cited as a solution with which customers have experienced significant success.<sup>3</sup> This widespread adoption by top-tier manufacturers underscores the reliability and effectiveness of HunterLab's technology in meeting the demanding color quality requirements of the potato chip market.

Beyond potato chips, HunterLab provides color measurement solutions for a diverse array of other snack foods, including cheese puffs, pretzels, cookies, and crackers.<sup>41</sup> These applications demonstrate the versatility of HunterLab spectrophotometers in addressing color quality across various food matrices that present different textural and structural characteristics. For instance, in the production of cheese puffs, precise color measurement is utilized as a key indicator of both taste and overall quality.<sup>41</sup> In the realm of cookies, smart spectrophotometers like the Aeros are capable of handling the often irregular textures of these baked goods, providing accurate and representative overall color measurements.<sup>43</sup> These examples illustrate that the principles and challenges of color control in non-uniform snack foods are common across different product categories, and HunterLab's technology is well-suited to address these needs.

The implementation of objective color measurement using HunterLab spectrophotometers in the snack food industry has led to significant benefits for

manufacturers. One key advantage is the reduction of waste. By enabling early detection and correction of color deviations during the production process, objective measurement helps prevent the manufacture of end products that do not meet the required color specifications.<sup>3</sup> Consistent color also directly contributes to an improvement in the overall quality of the product, leading to greater consumer satisfaction.<sup>9</sup> Perhaps most importantly for brand recognition and customer loyalty, maintaining precise color standards ensures brand consistency. Consumers come to expect a certain visual appearance from their favorite snack food brands, and consistent color reinforces this expectation, fostering trust and encouraging repeat purchases.<sup>8</sup> Therefore, investing in objective color measurement with instruments like the HunterLab Aeros offers tangible returns in terms of operational efficiency, enhanced product quality, and a stronger, more reliable brand presence in the marketplace.

## **IX. Conclusion: Achieving Superior Quality and Brand Recognition Through Precise Color Management**

In conclusion, the color of potato chips is a critical quality attribute that significantly influences consumer perception, brand loyalty, and overall product success. Implementing precise color measurement using a spectrophotometer like the HunterLab Aeros offers a powerful tool for ensuring both the quality and brand consistency of this popular snack food. The Aeros, with its non-contact measurement, automatic height positioning, and rotating sample platform, is specifically designed to address the challenges inherent in measuring the color of irregular, non-uniform samples like potato chips.

Achieving consistent and desirable color requires a comprehensive approach that extends beyond simply taking measurements. It necessitates a thorough understanding of the principles of spectrophotometry, the proper application of standard operating procedures for instrument use, and the careful selection and utilization of relevant color scales and indices such as CIE L\*a\*b\*, Hunter L, a, b, Agtron, and BCU. Furthermore, manufacturers must establish clear color ranges and target values based on consumer preferences and their unique brand specifications. Recognizing and controlling the influence of various processing parameters, including frying time and temperature, potato variety, maturity, storage conditions, and pre-treatments, is equally crucial.

By integrating data from the HunterLab Aeros into real-time monitoring systems and statistical process control, manufacturers can proactively manage color variations and ensure that their potato chips consistently meet the desired standards. The

experiences of numerous snack food companies, including leading potato chip producers who rely on HunterLab spectrophotometers, demonstrate the tangible benefits of this approach, including reduced waste, improved product quality, and enhanced brand recognition.

To optimize color management in potato chip production, it is recommended that manufacturers:

1. Establish clear and measurable color targets and tolerance limits using appropriate color scales based on consumer research and brand identity.
2. Implement a detailed standard operating procedure for the HunterLab Aeros, ensuring proper instrument setup, calibration, standardization, sample preparation, and data acquisition.
3. Integrate the Aeros into the production line for real-time color monitoring and feedback, allowing for timely adjustments to processing parameters.
4. Incorporate color measurement data into statistical process control systems to track process stability and identify sources of variation.
5. Work closely with potato suppliers to ensure a consistent supply of suitable potato varieties with low reducing sugar content and optimal storage conditions.
6. Optimize pre-treatment and frying processes to minimize color variability.
7. Regularly train personnel on the importance of color consistency and the proper use of the HunterLab Aeros.

By embracing a holistic strategy that combines advanced color measurement technology with a commitment to process control and quality at every stage of production, potato chip manufacturers can consistently deliver products that not only meet but exceed consumer expectations, thereby strengthening their brand reputation and ensuring long-term success in the marketplace.

### **Works cited**

1. Grading of Potato Chips According to Their Sensory Quality Determined by Color, accessed on March 31, 2025, [https://www.researchgate.net/publication/251264752\\_Grading\\_of\\_Potato\\_Chips\\_According\\_to\\_Their\\_Sensory\\_Quality\\_Determined\\_by\\_Color](https://www.researchgate.net/publication/251264752_Grading_of_Potato_Chips_According_to_Their_Sensory_Quality_Determined_by_Color)
2. Comprehensive Analysis of Physicochemical Properties and Sensory Attributes of Original-Cut Potato Chips in the Chinese Market - MDPI, accessed on March 31, 2025, <https://www.mdpi.com/2304-8158/13/24/4158>
3. What Is the Best Way to Measure the Color of Potato Chips? | HunterLab, accessed on March 31, 2025, <https://www.hunterlab.com/blog/the-best-way-to-measure-the-color-of-potato-chips/>

4. Impact of the Temperature Reconditioning of Cold-Stored Potatoes on the Color of Potato Chips and French Fries - PMC, accessed on March 31, 2025, <https://pmc.ncbi.nlm.nih.gov/articles/PMC10930719/>
5. Yield, Dry Matter, Specific Gravity and Color of Three Bangladeshi Local Potato Cultivars as Influenced by Stage of Maturity - Science Alert, accessed on March 31, 2025, <https://scialert.net/fulltext/?doi=jps.2015.108.115>
6. View of COLOUR AND CRISPNESS ASSESSMENT OF FORTY POTATO VARIETIES FOR PROCESSING INDUSTRY OF BANGLADESH, accessed on March 31, 2025, <https://epubs.icar.org.in/index.php/PotatoJ/article/view/59708/24609>
7. Reference chart for potato chip colour; colour score 1-6. - ResearchGate, accessed on March 31, 2025, [https://www.researchgate.net/figure/Reference-chart-for-potato-chip-colour-colour-score-1-6\\_fig1\\_321348534](https://www.researchgate.net/figure/Reference-chart-for-potato-chip-colour-colour-score-1-6_fig1_321348534)
8. Measure The Color Of Potato Chips With A Spectrophotometer, accessed on March 31, 2025, <https://www.hunterlab.com/en/industries/food/snack-foods/measuring-the-color-of-potato-chips/>
9. The Importance of Food Color Measurement | HunterLab, accessed on March 31, 2025, <https://www.hunterlab.com/blog/importance-of-food-color-measurement/>
10. United States Standards for Grades of Potatoes for Chipping - Agricultural Marketing Service, accessed on March 31, 2025, [https://www.ams.usda.gov/sites/default/files/media/Potatoes\\_for\\_Chipping\\_Standard%5B1%5D.pdf](https://www.ams.usda.gov/sites/default/files/media/Potatoes_for_Chipping_Standard%5B1%5D.pdf)
11. How Spectrophotometers Have Made the Visual Color Standards Reference Chart Obsolete for Potato Chips | HunterLab, accessed on March 31, 2025, <https://www.hunterlab.com/blog/how-spectrophotometers-have-made-the-visual-color-standards-reference-chart-obsolete-for-potato-chips/>
12. Spectrophotometers and the Color Standards Reference Chart for Potato Chip Analysis - MeasureColour, accessed on March 31, 2025, <https://measurecolour.com.my/wp-content/uploads/2019/03/Spectrophotometers-and-the-Color-Standards-Reference-Chart-for-Potato-Chip-Analysis.pdf>
13. Spectrophotometers and the Color Standards Reference Chart for Potato Chip Analysis, accessed on March 31, 2025, <https://www.hunterlab.com/blog/spectrophotometers-and-the-color-standards-reference-chart-for-potato-chip-analysis/>
14. www.hunterlab.com, accessed on March 31, 2025, [https://www.hunterlab.com/media/documents/120-Aeros\\_Solutions\\_Broch-Food-PotatoChips-04-22.pdf](https://www.hunterlab.com/media/documents/120-Aeros_Solutions_Broch-Food-PotatoChips-04-22.pdf)
15. Aeros | Non-Contact Spectrophotometer For Color Measurement - HunterLab, accessed on March 31, 2025, <https://www.hunterlab.com/en/products/benchttop-spectrophotometers/aeros/>
16. How to measure the color of Chips - YouTube, accessed on March 31, 2025, <https://www.youtube.com/watch?v=Xhh2CWYKrps>
17. www.hunterlab.com, accessed on March 31, 2025, [https://www.hunterlab.com/media/documents/Users\\_Manual\\_for\\_Aeros\\_1.1.pdf](https://www.hunterlab.com/media/documents/Users_Manual_for_Aeros_1.1.pdf)

18. How to Operate a Spectrophotometer - HunterLab, accessed on March 31, 2025, <https://staging.hunterlab.com/blog/how-to-operate-a-spectrophotometer/>
19. Standard Operating Procedures (SOP): Hunterlab Colourimeter - YouTube, accessed on March 31, 2025, <https://www.youtube.com/watch?v=c08qgdFG0-w>
20. COLOR MEASUREMENT OF FOOD PRODUCTS USING CIE L\*a\*b\* AND RGB COLOR SPACE - Journal of Hygienic Engineering and Design, accessed on March 31, 2025, <https://keypublishing.org/jhed/wp-content/uploads/2020/07/08.-Full-paper-Ivana-Markovic.pdf>
21. CIELAB color space - Wikipedia, accessed on March 31, 2025, [https://en.wikipedia.org/wiki/CIELAB\\_color\\_space](https://en.wikipedia.org/wiki/CIELAB_color_space)
22. Chips colour (L\*, a\*, b\*, chroma, hue angle) and crispness of forty potato varieties. - ResearchGate, accessed on March 31, 2025, [https://www.researchgate.net/figure/Chips-colour-L-a-b-chroma-hue-angle-and-crispness-of-forty-potato-varieties\\_tbl1\\_305220854](https://www.researchgate.net/figure/Chips-colour-L-a-b-chroma-hue-angle-and-crispness-of-forty-potato-varieties_tbl1_305220854)
23. Rapid Screen of the Color and Water Content of Fresh-Cut Potato ..., accessed on March 31, 2025, <https://pmc.ncbi.nlm.nih.gov/articles/PMC7022740/>
24. We Eat Snacks With Our Eyes First - Colors Matter: HunterLab Launches Snack Food Campaign, Providing State-of-the-Art Spectrophotometers to Ensure Vibrant Food Hues and Color Tones - PR Newswire, accessed on March 31, 2025, <https://www.prnewswire.com/news-releases/we-eat-snacks-with-our-eyes-first-colors-matter-hunterlab-launches-snack-food-campaign-providing-state-of-the-art-spectrophotometers-to-ensure-vibrant-food-hues-and-color-tones-301650373.html>
25. Hunter colour L (HC-L) and Yellow index (YO values corresponding to... - ResearchGate, accessed on March 31, 2025, [https://www.researchgate.net/figure/Hunter-colour-L-HC-L-and-Yellow-index-YO-values-corresponding-to-different-chip-colour\\_tbl1\\_321348534](https://www.researchgate.net/figure/Hunter-colour-L-HC-L-and-Yellow-index-YO-values-corresponding-to-different-chip-colour_tbl1_321348534)
26. Chip Potato Varieties for Commercial Production in Northeast Florida - UFDC Image Array 2, accessed on March 31, 2025, <http://ufdcimages.uflib.ufl.edu/IR/00/00/16/93/00001/CV28000.pdf>
27. M-SERIES III — Agtron Inc: Analyzers, Spectrophotometers, Colorimeters, accessed on March 31, 2025, <https://www.agtroninc.com/view-analyzers/m-series>
28. Products - ST Analytical, accessed on March 31, 2025, <https://www.agtron.eu/Products/>
29. A REFERENCE CHART FOR POTATO CHIP COLOUR FOR USE IN INDIA - ResearchGate, accessed on March 31, 2025, [https://www.researchgate.net/profile/Brajesh\\_Singh3/publication/321348534\\_A\\_REFERENCE\\_CHART\\_FOR\\_POTATO\\_CHIP\\_COLOUR\\_FOR\\_USE\\_IN\\_INDIA/links/5a1e3e2fa6fdccc6b7f879e4/A-REFERENCE-CHART-FOR-POTATO-CHIP-COLOUR-FOR-USE-IN-INDIA.pdf](https://www.researchgate.net/profile/Brajesh_Singh3/publication/321348534_A_REFERENCE_CHART_FOR_POTATO_CHIP_COLOUR_FOR_USE_IN_INDIA/links/5a1e3e2fa6fdccc6b7f879e4/A-REFERENCE-CHART-FOR-POTATO-CHIP-COLOUR-FOR-USE-IN-INDIA.pdf)
30. Potato Chip - Target, accessed on March 31, 2025, <https://www.target.com/s/potato+chip>
31. Potato Chip - Target, accessed on March 31, 2025, <https://www.target.com/c/chips-snacks-grocery/potato-chip/-/N-5xsy7Zyzj8n>

32. Effect of Slice Thickness and Frying Temperature on Color, Texture and Sensory Properties of Crisps made from Four Kenyan Potato Cultivars - Science Alert, accessed on March 31, 2025, <https://scialert.net/fulltext/?doi=ajft.2011.753.762>
33. www.researchgate.net, accessed on March 31, 2025, [https://www.researchgate.net/publication/275979519\\_Effect\\_of\\_Slice\\_Thickness\\_and\\_Frying\\_Temperature\\_on\\_Color\\_Texture\\_and\\_Sensory\\_Properties\\_of\\_Crisps\\_made\\_from\\_Four\\_Kenyan\\_Potato\\_Cultivars#:~:text=The%20effect%20of%20frying%20temperature,o%20C%20for%204%20mins.](https://www.researchgate.net/publication/275979519_Effect_of_Slice_Thickness_and_Frying_Temperature_on_Color_Texture_and_Sensory_Properties_of_Crisps_made_from_Four_Kenyan_Potato_Cultivars#:~:text=The%20effect%20of%20frying%20temperature,o%20C%20for%204%20mins.)
34. Potato chip color scores after frying tests in relation to different storage conditions., accessed on March 31, 2025, [https://www.researchgate.net/figure/Potato-chip-color-scores-after-frying-tests-in-relation-to-different-storage-conditions\\_tbl3\\_331780299](https://www.researchgate.net/figure/Potato-chip-color-scores-after-frying-tests-in-relation-to-different-storage-conditions_tbl3_331780299)
35. (PDF) Effect of Slice Thickness and Frying Temperature on Color, Texture and Sensory Properties of Crisps made from Four Kenyan Potato Cultivars - ResearchGate, accessed on March 31, 2025, [https://www.researchgate.net/publication/275979519\\_Effect\\_of\\_Slice\\_Thickness\\_and\\_Frying\\_Temperature\\_on\\_Color\\_Texture\\_and\\_Sensory\\_Properties\\_of\\_Crisps\\_made\\_from\\_Four\\_Kenyan\\_Potato\\_Cultivars](https://www.researchgate.net/publication/275979519_Effect_of_Slice_Thickness_and_Frying_Temperature_on_Color_Texture_and_Sensory_Properties_of_Crisps_made_from_Four_Kenyan_Potato_Cultivars)
36. Genotype × Storage Environment Interaction and Stability of Potato Chip Color, accessed on March 31, 2025, <https://horticulture.webhosting.cals.wisc.edu/wp-content/uploads/sites/20/2014/03/5-rak-20133.pdf>
37. Influence of potato composition on chip color quality | Semantic Scholar, accessed on March 31, 2025, <https://www.semanticscholar.org/paper/Influence-of-potato-composition-on-chip-color-Rodriguez-Saona-Wrolstad/942c83da43d5032db3d744572fef46ab7930102e>
38. Potato Variety Selection | College of Agricultural Sciences, accessed on March 31, 2025, <https://cropandsoil.oregonstate.edu/potatoes/potato-variety-selection>
39. Chip-Stock | Potato Goodness, accessed on March 31, 2025, <https://potatogoodness.com/wp-content/uploads/Chip-Stock-Broch-1.pdf>
40. HunterLab SpectraTrend HT – The Best Inline Color Measurement Solution for Baked Snack Food - YouTube, accessed on March 31, 2025, <https://www.youtube.com/watch?v=5sFWMQLcskg>
41. Measure The Color Of Cheese Puffs With A Spectrophotometer - HunterLab, accessed on March 31, 2025, <https://www.hunterlab.com/en/industries/food/snack-foods/measuring-the-color-of-cheese-puffs/>
42. Snack Foods Color Measurement Solutions - HunterLab, accessed on March 31, 2025, <https://www.hunterlab.com/en/industries/food/snack-foods/>
43. Smart Spectrophotometers Offer a Better Way to Measure the Color of Cookie Samples, accessed on March 31, 2025, <https://www.hunterlab.com/blog/smart-spectrophotometers-offer-a-better-way-to-measure-the-color-of-cookie-samples/>